

Introduction to CBC and Blood Smear Interpretation

Thanyaphong Na Nakorn, MD, PhD
Division of Hematology
Department of Medicine
Chulalongkorn University

1

Complete blood count

- The most common test used in clinical medicine
- Determine type and severity of blood cell abnormalities
- Nowadays, CBC is fully automated and highly reproducible.
- Correct interpretation of automated CBC can reduce rate of unnecessary blood smear examination
- Provide useful information for provisional diagnosis of RBC and WBC diseases

2

Component of automated CBC

- Blood count basic parameters: Hb, Hct, RBC, WBC, platelet.
- Red cell indices: MCV, MCH, MCHC, RDW
- WBC differentials
- Cytogram or Scattergram
- Reticulocyte count

3


Basic principles of automated blood cell analyzer

- Electrical impedance:
Coulter Counter , Cell-Dyn®
- Optical impedance and light scatter:
Technicon-H series
- VCS (volume, conductivity, light scatter) technology
Beckman-Coulter VCS, MAXM, STKS

4


Coulter Technology

- Cell count and size can be measured by electrical impedance



5

Hydrodynamic Focusing



The Bernoulli Effect




Viscous drag along walls.

Particles move to low pressure area



6

Technicon Technology



7

VCS Technology



8

Instrument	Methodology features
Coulter STKS	Impedance, conductivity, laser technology
Coulter GEN-S	Impedance, flow cytometry, conductivity, laser light scatter (reticulocyte analysis)
Sysmex SE-9000	Direct current (resistance to cell volume); radio-frequency (cell density, size)
Cobas Argos 5-diff (Roche)	Impedance, light absorption (halogen light source)
Technicon H-3 (Bayer)	Cytochemistry, flow technology (reticulocyte analysis)
Advia 120 (Bayer)	Cytochemistry, flow technology (reticulocyte analysis)
Cell Dyn 4000 (Abbott)	Impedance, laser light scatter


9

Red cell parameters

- Direct Measurement
 - Erythrocyte Concentration (RBC) $\times 10^6/\text{ml}$
 - Mean Corpuscular Volume (MCV) Femtolitre (fl)
 - Hemoglobin (Hb) Gram/decilitre (g/dl)
- Indirect Measurement
 - Hematocrit (Hct) = $\text{RBC} \times \text{MCV} / 10$ %
 - Mean Corpuscular Hemoglobin (MCH) = $\text{HB} \times 10 / \text{RBC}$ (pg)
 - Mean Corpuscular Hemoglobin Concentration (MCHC) = $\text{Hb} / \text{Hct} \times 100$ (g/dl)
 - Red Cell Distribution Width (RDW) %


10

RBC size and MCV




11

Red Cell Distribution Width



12

WBC: three-part differentials




Cell-Dyn 1700CS

13

WBC: Five-part differentials

- Neutrophil
- Lymphocyte
- Monocyte
- Eosinophil
- Basophil




Coulter VCS

14

WBC: Six-part differentials

- Neutrophil
- Lymphocyte
- Monocyte
- Eosinophil
- Basophil
- Large Unstained Cell (LUC)

Peroxidase channel




Technicon*H2

15

WBC: Six-part differentials

- Additional parameters in Technicon-H analyzer
 1. Mean peroxidase activity index (MPXI)
 2. Lobularity Index (LI)

Basophil/lobularity channel



Technicon*H2

16



Platelet parameters

- Platelet count
- Mean Platelet Volume (MPV)
- Platelet Distribution Width (PDW)
- Plateletcrit (Pct)

17


Reticulocyte count

- Reticulocyte = non-nucleated RBC with polyribosomal RNA as stained by supravital stain (new methylene blue or brilliant cresyl blue)
- Polychromasia underestimates reticulocytes
- Three methods of reticulocyte enumeration
 - Manual count on slide per 1,000 RBC
 - Automated CBC with reticulocyte counter (Coulter VCS, Cell-Dyne 4000, Technicon-H3)
 - Flow cytometry with fluorescent dyes




18

Report from Coulter™




19

Report from Technicon-H™



20



21


Clinical Utility of automated CBC

- Differential diagnosis of RBC disorders
 - Microcytic RBC
- Use WBC differential cytogram to guide WBC morphology and phenotypes

22

What is your diagnosis?

WBC 5.00 $\times 10^3/\mu\text{L}$
 RBC 3.56 $\times 10^6/\mu\text{L}$
 Hb 6.0 g/dL
 HCT 21.0 %
 MCV 59.0 fL
 MCH 17.0 pg
 MCHC 28.8 g/dL
 RDW 19.4 %
 HDW 4.02 g/dL
 PLT 427 $\times 10^3/\mu\text{L}$



23

Differential diagnosis of anemia using MCV and RDW

	Low MCV	Normal MCV	High MCV
RDW <15	Thalassemia trait Heterozygous HbE, HbC, etc. Anemia of chronic disease (ACD)	ACD Heterozygous HbS, HbCS, HbE, etc. Hereditary spherocytosis Acute hemorrhage	Aplastic anemia MDS Myeloma Liver disease Hyperthyroidism
RDW >15	Iron deficiency anemia Thalassemia intermedia Sideroblastic anemia Severe ACD RBC fragmentation	Early or combined nutritional deficiency Myelodysplasia Myelophthisis Sickle cell anemia or Homozygous HbCS	B12 deficiency Folate deficiency AIHA Drugs: HU, ARV, AZA, etc.

24

Factors known to cause spurious laboratory results in hematology analyzers.

Parameter	Spuriously increased	Spuriously decreased
RBCs	WBC >50 000/mm ³	Clotting
Hb	hyperlipidemia, hyperbilirubinemia	Clotting
MCV	Cold agglutinins, hyperglycemia, WBC >50 000/mm ³	Cryoglobulins
MCHC	Hyperlipidemia, cold agglutinins	WBC >50 000/mm ³
RDW	Post transfusion	
WBCs	Nucleated red cells, platelet clumps, unlysed red cells, cryoglobulins	Clotting
Platelets	WBC fragmentation, severe microcytosis, cryoglobulins	Satellitism, clumping

25

Basic Principles for Blood Smear Interpretation

- Assess quality of smears
 - Specimen preparation & staining
- Estimate cell numbers
 - RBC: evenly dispersed with minimal intercellular space
 - WBC: 10-20/LPF
 - Platelets: 7-20/OF
- Determine predominant cell populations
- Carefully examine cellular morphology


26

RBC disorders

- Hypochromic microcytic anemia
 - Iron deficiency anemia
 - Thalassemia and hemoglobinopathy
- Macrocytic anemia
 - Megaloblastic anemia
 - Non-megaloblastic macrocytic anemia
- Hemolytic anemia
 - Immune hemolytic anemia: AIHA, DHTR
 - Microangiopathic hemolytic anemia (MAHA)
 - Red cell enzymopathies: G-6-PD deficiency
 - RBC membrane defects: spherocytosis, ovalocytosis, elliptocytosis, stomatocytosis
- RBC inclusion bodies and parasites


27

Area of blood film examination




28


RBC Size and MCV




29




30




31




32










35




36







41




42




43



44




45




46




47




48




49




50




51




52




53




54



55




56




57




58




59




60




61




62



63



64




65



WBC disorders

- Leukopenia
 - with absolute neutropenia: bone marrow failure, agranulocytosis
 - with atypical lymphocytes: viral infection, chronic lymphoproliferative disorders
 - with immature myeloid cells: acute leukemia, MDS or myelophthisis
- Leukocytosis
 - Reactive leukocytosis: leukemoid reaction
 - Acute leukemia: AML vs. ALL
 - Chronic myeloproliferative disorders
 - Chronic lymphoproliferative disorders
- Leukoerythroblastosis


67



68




69




70




71





72




73




74







77





78



79




80




Platelet disorders

- Quantitative disorders
 - Isolated thrombocytopenia: Immune vs. non-immune
 - Thrombocytopenia associated with other hematologic abnormalities
 - Thrombocytosis
- Qualitative disorders
 - Giant platelets (megathrombocytes)
 - Platelet inclusion or granule abnormality
 - Bizarre in shape and size
 - Megakaryocytes or megakaryoblasts


83




84




85




86



87




88




89




90



91



92



93